# Momini

## NebuBunny

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Create happy beginnings

### ● Please read this guidebook carefully before this unit

Your Nebulizer is intended for use in the treatment of asthma, COPD and other respiratory ailments in which an aerosolized medication is required during therapy. Please consult with your physician and/or pharmacist to determine if your prescription medication is approved for use with this Nebulizer. This Nebulizer is intended for single patient use.


### ● 01. INTRODUCTION

Thank you for purchasing the Compressor Nebulizer. It is a compact medical device designed to efficiently deliver physician prescribed medication to the bronchial lung passages. With proper care and use, it will provide you with many years of reliable treatment.

This product is developed for the successful treatment of asthma, allergies and other respiratory disorders. It creates a stream of air that travels through clear tube to the nebulizer. When air enters the nebulizer, it will convert the prescribed medication into aerosol mist for easy inhalation.

Your Compressor Nebulizer should be used under the supervision of a licensed physician and/or a respiratory therapist. We encourage you to thoroughly read this guidebook to learn about the features of this product. Any use of this product other than its intended use should always be avoided.

### ● 02. PRODUCT IDENTIFICATION

### ● 03. IMPORTANT SAFEGUARDS

**Note:** Read all instructions carefully before use. The following basic precautions are needed when using an electrical product:

**Caution:** Failure to read and observe all precautions could result in personal injury or equipment damage.

### Product cautions:

1. To avoid electrical shock: Keep the unit away from water • Do not immerse the power cord or the unit in liquid • Do not use while bathing • Do not reach for a unit that has fallen into water – immediately unplug the unit.
2. Never operate the unit if it has any damaged parts (including power cord), if it has been dropped or submerged in water. Promptly send it to a service center for examination and repair.
3. The unit should not be used where flammable gas, oxygen or aerosol spray products are being used.
4. Disconnect the unit from the electrical outlet before cleaning, filling and after each use.
5. Do not use other attachments unless recommended by manufacturer.

### Operating Cautions:

1. Connect this product to an appropriate voltage outlet for your model.
2. Do not run this product unattended.
3. Never operate if this unit has a damaged cord or plug, if it has been dropped, or in any way, submerged into water, if it does not work properly. Return it to a service center for repair.
4. If any abnormality occurs, discontinue use immediately until the unit has been examined and repaired.
5. Always unplug the product immediately after use.
6. Never block the air openings of the main unit or place it where the air openings may be obstructed.

### Storage Cautions:

1. Do not store the unit under direct sunlight, high temperature or humidity.
2. Keep the unit out of reach of small children.
3. Always keep the unit unplugged while not in use.


### Cleaning Cautions:

1. Do not immerse the unit in water. It may damage the unit.
2. Disconnect the unit from the electrical outlet before cleaning.
3. Clean all necessary parts after each use as instructed in this guidebook.

### ● 04. OPERATING YOUR COMPRESSOR NEBULIZER

**Note:** Prior to initial operation, the nebulizer should be thoroughly cleaned.

1. Place your Compressor Nebulizer on a flat and stable surface. Be sure that you can easily reach the controls when you are seated.
2. Gently twist the top part of the nebulizer counter-clockwise to disassemble the nebulizer.
3. Fill the bottom section of the nebulizer with the medication prescribed by your physician. Be sure that the cone is put inside the bottom section.
4. Gently twist the top part clockwise to reassemble the nebulizer. Be

sure the two sections fit well.

5. Attach one end of air tube to the base of the nebulizer.
  6. Attach the other end of air tube to the tube connector located on the front of the compressor.
  7. Attach mouthpiece or mask of your choice to the top section of nebulizer.
  8. Plug the power cord into an appropriate electrical outlet. Make sure at this stage, the power switch is at "Off" status.
  9. Press power switch to begin your prescribed treatment.
- Important:** The compressor motor has a thermal protector, which will shut off the unit before the unit is overheated. When the thermal protector shuts the unit off, please:
- a. Switch off the unit.
  - b. Unplug the unit from the electrical outlet.
  - c. Wait 30 minutes for the motor to cool down before another treatment. Make sure the air openings are not obstructed.
10. When treatment is finished, shut off the unit and unplug it from the electrical outlet.

### ● 05. CLEANING

Cleaning & Disinfection for the Nebulizer's Accessories (including: nebulizer cup, mouthpiece, child mask, adult mask, nose piece, air tube)

#### After Use the Nebulizer Accessories

1. Clean the Nebulizer's Accessories with distilled water.
2. Completely immerse the Nebulizer's Accessories into approx 75% medicinal alcohol.
3. Soak for approx 30 minutes.
4. Wash and rinse the Nebulizer's Accessories with about 50 warm distilled water.
5. Shake out water dry them with a clean, soft towel or leave them air-dry in a clean place.
6. Make sure that the Nebulizer's Accessories and other cleaned parts are completely dry then put all these parts in a dry sealed container.

#### Cleaning the compressor

1. Wipe daily with a damp cloth.

2. Do not use any powdered cleaners or soap pads, which may damage the finish.

#### Filter change

1. Do not use cotton or any other materials. Do not wash or clean the filter. Only use filters supplied by your distributor. And do not operate without a filter.
2. Change the filter every 30 days or when the filter turns gray.
3. Changing procedure:
  - A. Remove the filter cover.
  - B. Replace the used filter with a new one.
  - C. Replace the filter cover.

#### 06. SPECIFICATIONS

|                                |                                                                                         |
|--------------------------------|-----------------------------------------------------------------------------------------|
| Model No.                      | MCN-S600MA                                                                              |
| Power                          | 230V/50HZ                                                                               |
| Power Consumption:             | 180VA                                                                                   |
| Nebulizer medication Capacity: | 7ml                                                                                     |
| Partical Size:                 | Approx.(78.25±13.08)% within 0.5 to 5um range                                           |
| MMAD:                          | 2.52±0.89um                                                                             |
| Average Nebulization Rate:     | ≥0.2ml/min                                                                              |
| Noise Level:                   | Below 65dBA (DIN EN 13544-1)                                                            |
| Compressor Pressure Range:     | 30 ~ 50 Psi / 206Kpa ~ 345Kpa                                                           |
| Operation Pressure Range:      | 9.5 ~ 16 Psi / 65.5Kpa ~ 110Kpa                                                         |
| Liter Flow Range:              | 6 ~ 10 lpm                                                                              |
| Operation Temperature Range:   | 10°C ~ 40°C (50°F ~ 104°F)                                                              |
| Operation Humidity Range:      | 10% ~ 95% RH                                                                            |
| Storage Temperature Range:     | -20°C ~ 70°C (-4°F ~ 158°F)                                                             |
| Storage Humidity Range:        | 10% ~ 95% RH                                                                            |
| Atmospheric pressure:          | 500 ~ 1060hPa                                                                           |
| Dimension (L x W x H)          | 143*136*200mm                                                                           |
| Weight                         | 1.34 kg                                                                                 |
| Standard Accessories           | Nebulizer Cup, Air Tube, Adult Mask, Child Mask, Mouthpiece, Nose Piece, Filters (5PCS) |

#### Operation Mode

Protection against electric shock:

Continuous

Class II Equipment

Type BF applied parts:

Mouthpiece, mask

CE 1639

#### 07. EU DECLARATION OF CONFORMITY

Dongguan Aidsy Machinery & Electronic Equipment Co., Ltd.

File No.: AIDISY-DCN-CE-02  
Version: V5.0

#### EC Declaration of Conformity

| Manufacturer:                                                                                                                                  | Whose Single Authorized Representative:           |
|------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------|
| Dongguan Aidsy Machinery & Electronic Equipment Co., Ltd.                                                                                      | MedPath GmbH                                      |
| Address:                                                                                                                                       | Address:                                          |
| Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhouwu, Dongcheng District Dongguan City, Guangdong Province, 523118 P.R. China | Mies-van-de-Rohe-Strasse 8, 80807 Munich, Germany |
| SRN: CN-MF-000026848                                                                                                                           | SRN: DE-AR-000000087                              |

We, the manufacturer, declare on our sole responsibility that following products

| Product name                               | Model                                                                                                                                                                                                                                               | Basic UDI-DI                       | Intended purpose                                                                                                                                                          |
|--------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Compressor Nebulizer for Conscious Patient | MCN-D5600A, MCN-D5600B, MCN-D5600C, MCN-D5600D, MCN-D5600E, MCN-S600A, MCN-S600B, MCN-S600C, MCN-S600D, MCN-S600E, MCN-S600F, MCN-S600G, MCN-S600MA, MCN-S600MB, MCN-S600MC, MCN-S600MD, MCN-S600ME, MCN-S600MF, MCN-S600MG, MCN-S600MH, MCN-S600MI | Basic UDI-DI: 69571303/MCN 0001 GX | The Compressor Nebulizer is intended to provide air to the pneumatic nebulizer in order to aerosolize medications for inhalation by the patient for respiratory disorders |

meet the provisions of Regulations (EU) 2017/745.

The medical device has been assigned to class IIa according to rule 20 in Annex VIII of MDR 2017/745. It bears the mark

CE 1639

The product concerned has been designed and manufactured under a quality management system according to Annex IX of Regulations (EU) 2017/745.

Compliance of the designated product with the Annex IX of Regulation (EU) 2017/745 has been assessed and certified by the Notified Body

SGS Belgium NV  
SGSHouse Noorderlaan  
87 2030 Antwerp Belgium

|                  |                                   |
|------------------|-----------------------------------|
| Certificate No.: | Certificate No.: CN1941056        |
| Issue date:      | SGS Belgium NV Notified Body 1639 |
| Expiry date:     | 24 MAY 2024                       |

The above mentioned declaration of conformity is exclusively under the responsibility of

|          |                                                                                                                                                |
|----------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Company: | Dongguan Aidsy Machinery & Electronic Equipment Co., Ltd                                                                                       |
| Address: | Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhouwu, Dongcheng District Dongguan City, Guangdong Province, 523118 P.R. China |

2024/05/01

Place, date

Legally binding signature, Function

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**[EN] Rules for handling electrical and electronic equipment:** Waste electrical and electronic equipment as well as batteries and accumulators cannot be mixed with other waste. This is indicated by the crossed out bin symbol on the product packaging. Electrical and electronic equipment as well as batteries and accumulators may contain hazardous substances, mixtures and components harmful to the environment and human health. Improper handling or damage may result in damage to health or environmental pollution. Proper household handling of waste electrical and electronic equipment and batteries and accumulators positively contributes to the reuse and recovery, including recycling, of waste equipment. Waste electrical and electronic equipment as well as batteries and accumulators can be returned to local selective waste collection points. The distributor is obliged to collect used equipment from households free of charge at the point of sale, provided that the used equipment is of the same type and performed the same functions as the equipment sold. When delivering equipment intended for households to the buyer, the distributor is obliged to collect used equipment

from households free of charge at the place of delivery of the equipment, provided that the used equipment is of the same type and performed the same functions as the delivered equipment. A distributor running a retail unit with a sales area within the meaning of Art. 2 point 19 of the Act of 27 March 2003 on spatial planning and development of at least 400 m<sup>2</sup> devoted to the sale of equipment intended for households, is obliged to accept waste equipment from households, free of charge, in this or in its immediate vicinity, none of the external dimensions exceeds 25 cm, without the need to purchase new household equipment. The product was introduced to the market on August 13, 2005 r.

**Manufactured by:** Dongguan Aidsy Machinery & Electronic Equipment Co., Ltd., Part B, 3rd F, Block A, Wentang Industrial Park, Longhua Road, Zhoulu, Dongcheng District, Dongguan City, Guangdong Province, 523118, P.R. China

**Authorized representative:** MedPath GmbH, Mies-van-der-Rohe-Strasse 8, 80807 Munich, Germany

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**EN |** More information about the product is available on the manufacturer's website at:

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