# TAKE CARE!

babyon

|           |   |     |
|-----------|---|-----|
| <b>GB</b> | Touch-free electronic thermometer         | 2   |
| <b>PL</b> | Termometr elektroniczny bezdotykowy       | 12  |
| <b>BG</b> | Електронен безконтактен термометър        | 22  |
| <b>CZ</b> | Elektronický bezkontaktní teploměr        | 32  |
| <b>DE</b> | Kontaktloser elektronischer Thermometer   | 42  |
| <b>ES</b> | Termómetro electrónico sin contacto       | 52  |
| <b>FR</b> | Thermomètre électronique sans contact     | 62  |
| <b>GE</b> | უკონტაქტო ელექტრო თერმომეტრი              | 72  |
| <b>GR</b> | Ανέπαφο ηλεκτρονικό θερμόμετρο            | 82  |
| <b>HU</b> | Érintés nélküli elektromos hőmérő         | 92  |
| <b>LT</b> | Elektroninis bekontaktnis termometras     | 102 |
| <b>LV</b> | Elektroniskais bezkontakta termometrs     | 112 |
| <b>NL</b> | Contactloze elektronische thermometer     | 122 |
| <b>PT</b> | Termómetro eletrónico sem contacto        | 132 |
| <b>RO</b> | Termometru digital non-contact            | 142 |
| <b>RU</b> | Термометр электронный бесконтактный       | 152 |
| <b>SE</b> | Elektronisk beröringsfri termometer       | 162 |
| <b>SK</b> | Elektronický bezkontaktný teplomer        | 172 |
| <b>TR</b> | Elektronik temasız/kızılıtesli termometre | 182 |
| <b>UA</b> | Термометр електронний безконтактний       | 192 |
| <b>XS</b> | Elektroniski bezkontaktni termometr       | 202 |

# Non-contact Infrared Thermometer

The non contact Infrared body thermometer operating Instructions intend to provide the necessary information for proper operation of HTD8813 thermometer model. Only body mode was reviewed and certified by notified body. General knowledge of Infrared thermometer and an understanding of the features and functions of the HTD8813 thermometer model are prerequisites for proper use. The non contact infrared body thermometer is a medical device, and can be used repeatedly, which using life is 5 years. Please read the manual first before using it, if not fully understand the usages, please stop using the thermometer. Do not operate any of the models HTD8813 thermometer without completely reading and understanding these instructions.

## 1. GENERAL INFORMATION

### INTENDED USE

Babyono non-contact infra-red thermometers were designed to check and monitor human body temperature by users themselves, in domestic setting, and primary care setting, for screening purposes. To learn how to operate the thermometer, please read the operating manual before using it.


### DESCRIPTION OF THE NON-CONTACT INFRA-RED THERMOMETER

The Babyono non-contact Infra-red thermometer is a manual, battery-powered device for measuring human body temperature (on the forehead). The device uses an infra-red sensor to measure the temperature. The IR sensor generates various signals when measuring temperatures of different objects or when placed in ambient with varying temperature. The ASIC transforms the infra-red signal into a digital value and shows on the device's LCD display.

### DESCRIPTION OF CONTROL ELEMENTS, INDICATORS AND SYMBOLS

Fig. 1: Infra-red thermometer, part no. 613

1. LCD display
2. Battery cover
3. On / Off
4. SET button
5. MEMO button
6. MODE button
7. IR sensor
8. Surface temperature measurement mode
9. Body temperature measurement mode
10. Readout
11. Result indicator
12. Low battery indicator
13. Sound on / off indicator
14. Saved readout number
15. Fahrenheit degrees
16. Celsius degrees
17. Memory indicator

**TECHNICAL SPECIFICATION**

|  |  |
|--|--|
| Measuring unit                           | °C / °F  |
| Operating mode                           | Adjusted mode (body temperature)   |
| Direct mode (surface temperature)        |  |
| Measuring site                           | Forehead   |
| Additional suggested measuring point     | behind the ear   |
| Reference Body Site                      | Axillary   |
| Rated output range                       | Body mode: 34.0 °C ~ 43.0 °C / 93.2 °F ~ 109.4 °F<br>Surface mode: 0 °C ~ 100 °C / 32 °F ~ 212 °F  |
| Range                                    | Body temperature measurement mode: 34.0 °C ~ 43.0 °C / 93.2 °F - 109.4 °F<br>Surface temperature measurement mode: 0 °C ~ 100.0 °C / 32 °F - 212 °F  |
| Accuracy                                 | Body temperature measurement mode:<br>34.0 °C ~ 34.9 °C ± 0.3 °C / 93.2 °F - 94.8 °F ± 0.5 °F<br>35.0 °C ~ 42.0 °C ± 0.2 °C / 95.0 °F - 107.6 °F ± 0.4 °F<br>42.1 °C ~ 43.0 °C ± 0.3 °C / 107.8 °F - 109.4 °F ± 0.5 °F<br>Surface temperature measurement mode:<br>± 2 °C / ± 3.6 °F   |
| Display resolution                       | 0,1 °C / 0,1 °F  |
| Three-colour backlight<br>(colour codes) | 35.5 °C - 37.3 °C / 95.9 °F - 99,1 °F Green (normal temperature);<br>37.4 °C - 38.0 °C (alarm point) / 99.3 °F - 100.4 °F Yellow (light fever)<br>38.1 °C - 43.0 °C / 100.6 °F - 109.4 °F Red (high fever)<br>Caution:<br>1. The colour codes only work in body temperature measuring mode.<br>2. The backlight is always green when measuring the temperature of surfaces.<br>3. In body mode 34.0-35.4°C is with green backlight.<br>Automatic power off time ≤ 18 s<br>Measuring time ≤ 2 s<br>Measuring distance 1 cm - 5 cm |
| Memory slots                             | 50   |
| Power supply requirements                |  |
| Batteries                                | 1.5 V (AAA) alkaline x 2 (IEC LR03)  |
| Operating voltage range                  | 2.6 V ~ 3.6 V  |
| Ambient conditions                       |  |
| Operating conditions                     | Operating temperature: 15 °C ~ 40 °C (59 °F ~ 104 °F),<br>relative humidity ≤ 85%, atmospheric pressure: 70 kPa - 106 kPa  |
| Transport and storage conditions         | Storage temperature: - 20 °C to + 55 °C / - 4 °F to + 131 °F,<br>relative humidity ≤ 93%, atmospheric pressure: 70 kPa - 106 Pa  |
| Dimensions and weight                    |  |
| Weight (w/o batteries)                   | 90 g   |
| Dimensions                               | L: 138 mm X W: 95 mm X H: 40 mm  |
| Compliance                               |  |
| Item                                     | Compliance with  |
| Device classification                    | Safety standard EN 60601-1: 2006+A1:2013, EN 60601-1-2: 2015   |
| Protection type                          | Internally (battery) powered equipment   |

## CALCULATED VALUES OF THE INDICATORS ACCORDING TO ISO 80601-2-56

| Indicators                         | Calculated value |
|------------------------------------|------------------|
| Clinical bias, $\Delta_{cb}$       | -0.027           |
| Standard deviation, $\sigma_j$     | 0.14             |
| Limits of agreement, LA            | 0.26             |
| Clinical repeatability, $\sigma_r$ | 0.07             |

## 2 - OPERATION

### 2.1 BATTERY INSTALLATION

Caution: A non-contact infra-red thermometer with discharged batteries will not operate. New batteries must be installed in such a case.

- 1) Pull down and open the battery cover in the thermometer's bottom part.
- 2) Place two AAA batteries inside, observing the "+" and "-" polarity marks.
- 3) Close the battery cover.

### 2.2 OPERATING METHOD

Please read the warnings in the manual before operating the device.

- To measure human body temperature, aim the thermometer at the middle of the forehead, maintain the distance and press the On/Off button - the measurement will take place immediately.
- The ambient temperature where the device is being used should be stable - do not approach large ventilators, AC vents, etc.
- When moving from a place with lower ambient temperature to a place with a higher temperature, remain for at least 5 minutes in the new place to allow the temperatures to balance out.
- Wait at least 1 second between measurements. When taking 5 (or more) consecutive measurements, wait at least 30 seconds between such series.
- The green backlight means the device is ready to use. The yellow backlit range of temperatures from 37.4 °C to 38.0 °C means light fever warning. Please take relevant steps if this body temperature is in this range. The red backlit temperature of over 38.1 °C means a fever. Please take relevant steps to lower body temperature or consult a doctor.
- Do not use the thermometer in places with excessive sunlight.
- If the measurement taken on the forehead is too low, measure the temperature behind the hair.

### 2.3 GENERAL CONFIGURATION AND OPERATION

- Starting the measurement

1. Use the On / Off button to activate the thermometer. After activation, the thermometer will take 2 seconds to perform a self-test.
2. Aim the thermometer at the forehead of the examined person, maintain the distance, press the On / Off button to measure the temperature and read the value on the LCD display.

Caution: 1) Once the display is off, the device makes a beep to signal the end of the measurement. At the same time, the LCD display will show the measured value and the background will highlight in either red, yellow or green.

- 2) To ensure accurate measurements, wait 30 seconds between making sets of five consecutive measurements.

- Selecting the operating mode

After activating the device, press the MODE button to select the body temperature measurement mode or the surface temperature measurement mode.

The body temperature measurement mode is for measuring the temperature of human body and the surface temperature measurement mode serves to measure the temperature of objects. (By default, the device operates in the body temperature measurement mode).

- Restoring and deleting the saved data

The 50-cell memory stores the result of the last measurement taken before deactivating the thermometer.

- 1) To display the stored measurements, briefly press the MEMO button when the device is being activated or deactivated.
- 2) An empty memory cell shows “---°C” or “---°F”.
- 3) The measured results can be stored in the device’s memory. The memory can store 50 results and once all cells are used, the oldest result is replaced with the latest result.
- 4) To fully wipe the memory, press the MEMO button when the device is being activated and keep it depressed, until “CLR” is shown on the LCD display.

## 2.4 CONFIGURING THE PARAMETERS

To satisfy the needs of various groups of patients or individual patients, the device’s colour codes can be customized.

To modify measurement parameters, press and keep depressed the SET button.

### 1) Selecting the measuring units - F1

When the device is being activated, press and keep depressed the SET button, then select F1, press and keep depressed the MODE button to switch between Celsius and Fahrenheit degrees, and then press the SET button to confirm the selection (Celsius degrees are selected by default).

### 2) Light fever alarm - F2

After selecting F1, press the SET button to go to F2, then press the MODE button to decrease the value by 0.1 °C, or press the MEMO button to increase the value by 0.1 °C; to adjust the value more quickly, press and keep the relevant button depressed, then press the SET button to save the selected value. (The default value is 38.1 °C)

### 3) Changing sound signals - F3

After selecting F2, briefly press the SET button to go to F3, then press the MODE button or the MEMO button to change the sound signal, and then press the SET button to confirm the new setting. (By default, there is a sound signal when the device is being activated).

### 4) General temperature offset value - F4

This value serves to adjust the range of temperature, depending on the patient’s unique characteristics or ambient conditions significantly affecting the measured temperature.

After selecting F3, briefly press the SET button to go to F4, then press the MEMO button to increase the value by 0.1 °C, or press the MODE button to decrease the value by 0.1 °C; to adjust the value more quickly, press and keep the relevant button depressed, then press the SET button to save the selected value.

Parameter adjustment range: from - 5 °C to + 5 °C (the default value is 0).

### 5) Exiting the configuration mode

After selecting F4, press the SET button to deactivate the LCD display and leave the settings menu.

- Activating and deactivating the sound signal

When the device is being activated, briefly press the SET button to turn the sound signal on or off.

- Restoring default settings

When the device is being activated, press and keep depressed the MODE button, until “rst” is shown on the LCD display. After 2 seconds, the default settings for positions F1 - F4 will be restored.

### 3 - TROUBLESHOOTING

| MESSAGE  | SITUATION   | SOLUTION  |
|--|---|---|
|  <p>Body <math>\text{^{\circ}C}</math> Body <math>\text{^{\circ}C}</math></p> | <p>The measured temperature is outside the typical range of human body temperature.<br/>(34,0 ~ 43,0 <math>\text{^{\circ}C}</math> / 93,2 <math>\text{^{\circ}F}</math> ~ 109,4 <math>\text{^{\circ}F}</math>).</p> | <p>Make sure the measurement is taken on the forehead, not on other body parts.</p>   |
|  | <p>The measurement should be taken from a distance of 1 - 5 cm.</p>   | <p>The optimal measuring distance is 1 cm.</p>  |
|  | <p>The examined person's forehead is covered with hair, sweat, anti-fever pads, etc.</p>  | <p>The examined person should wait, in leisure, 15 minutes, before having their temperature measured.</p>   |
|  | <p>Some individuals may have higher natural body temperature than the rest of the population.</p>   | <p>High temperature</p>   |
|  | <p>F4 improperly set general temperature offset</p>   | <p>Adjust the temperature offset value</p>  |
|  <p>Body <math>\text{^{\circ}C}</math></p>                                    | <p>The operating temperature exceeds the specified limits.</p>  | <p>Go to a room where the temperature is within the specified limits and wait 30 minutes before taking measurements.</p>  |
|   | <p>The LCD display is blinking and the device deactivates automatically.</p>  | <p>Replace the batteries. Have the device repaired.</p>   |
|   | <p>Discharged batteries. Temperature cannot be measured.</p>  | <p>Install new batteries.</p>   |
|   | <p>Excessive ambient temperature changes.</p>   | <p>Wait until the ambient temperature is stabilized.</p>  |
|   | <p>(1) Power supply deactivated.<br/>(2) Improperly installed batteries.<br/>(3) Discharged batteries.<br/>(4) No information shown on the LCD display.</p>   | <p>(1) Press the On / Off button again.<br/>(2) Check battery polarization.<br/>(3) Replace the batteries.<br/>(4) Contact the manufacturer of the service.</p> |

### 4 - BATTERY REPLACEMENT

1. Open and remove the battery cover, as shown on the cover. Make sure the device is off, before replacing the batteries.
  2. Remove the discharged batteries and replace with two new AAA batteries. Make sure the batteries are inserted as shown on the cover.
  4. Slide the battery cover back into place.
- Do not throw away the used batteries together with domestic waste. Old batteries should be disposed of in specially designated collection places.
5. If the device is still not working after replacing the batteries, make sure that the batteries have been inserted properly. Remove the barriers, wait 30 seconds and reinstall them.

### WARNING

The batteries must not be charged, disassembled or placed in fire.

1. The typical life of new, fresh batteries is 2,000 measurements, each lasting 18 seconds.

2. Only the recommended batteries should be used. Regular batteries should not be recharged. No batteries should be placed in fire.

3. Remove the batteries from the device if you intend to store it for an extended time.

## 5 - CLEANING, MAINTENANCE AND STORAGE

The sensor's lens is very delicate.

The lens must be protected from damage and dirt.

A clean, dry piece of cloth should be used to clean the device and its LCD display. The device should not be cleaned with solvents or submerged in water or other liquids.

The thermometer should be always stored in places where the ambient temperature and humidity are as specified in Item 1.

The thermometer should be stored in a dry, dust-free place. The thermometer should not be exposed to direct sunlight, high temperature or humidity, because such conditions may adversely affect its operation. If the ambient temperature changes drastically, e.g. when moving the thermometer from a cold place to a hot place, or if the ambient temperature is from 15 °C to 40 °C, wait for 30 minutes before using the device.

## 6 - DISPOSAL

1) Do not throw away the used batteries together with domestic waste. Old batteries should be disposed of in specially designated collection places.

2) Do not throw away the used device together with domestic waste. Please find out where this device can be safely disposed of. The local legal regulations currently in force should be observed.

## 7 - WARRANTY

The non-contact infra-red thermometer is covered with the manufacturer's 2-year warranty, starting from the date of purchase.

The warranty is void in the following situations:

- The serial number has been removed from the device or becomes illegible.
- The device was damaged due to being improperly connected to other devices.
- The device was damaged in an accident.
- The device was modified by the user without the manufacturer's previous consent.
- The batteries and the packaging are not covered by the warranty.
- When filing a warranty claim, the user shall present a warranty card indicating the date of purchase and the seller's stamp (including their name and address). Please request to stamp your warranty card when purchasing the device. The defective product should be delivered to a service. Defective products not covered by warranty can be repaired at a fee.

### CAUTION:

1. If there are any problems with the device, e.g. when setting it up, maintaining or operating, please contact Babyono's service department. The device should not be opened or repaired by the user.
2. All cases of non-standard operation should be submitted to Babyono's service department.
3. The device should be calibrated every two years or when subject to shock.
4. The device was designed for use by patients. The patient may make measurements and replace batteries under normal operating conditions and may also maintain the device as indicated in this manual.

### WARNING

- A defective or improperly operating thermometer should not be used for any purpose.
- If the ambient temperature changes drastically, e.g. when moving the thermometer from a cold place to a hot place, or if the ambient temperature is from 15 °C to 40 °C, wait for 30 minutes before using the device.
- Remove the batteries from the thermometer if it is not going to be used for a longer time.
- The device is not waterproof and should not be submerged in water or other liquids. Follow the instructions under Item 5 when cleaning or disinfecting the device.
- The IR sensor should not be touched.
- The thermometer should not be used if anti-fever pads or other methods of reducing body temperature are used.
- The thermometer should not be moved close to strong electrostatic or magnetic field, because this may affect the accuracy of measurements.
- To avoid damaging the device, old and new batteries should not be used simultaneously.
- If the examined person's forehead is covered with hair, sweat or a headscarf, the measurement may be inaccurate.
- The results of measurements made with the device are only for reference purposes. If in doubt, other methods of measuring temperature should be used.
- The device should be stored out of reach of children and pets. If the device is not going to be used, it should be stored in a dry room, protected against excessive humidity, heat, dust and direct sunlight. No heavy items should be placed on the device.
- The device should not be dropped, disassembled or modified.
- The device should not be used if it is damaged or behaves abnormally.
- The device contains sensitive components and should be handled with care. The storage and operating conditions described in Item 1 "Technical specification" should be observed.
- No potential allergic reactions are known as per the current state of technical knowledge.
- This thermometer will irreplaceable the diagnostic in hospitals.



For further information contact:

Manufacturer: Hetaida Technology Co., Ltd. Add: 4F, BaiShiDa High-Tech Park, XiangDong Industrial Area, DaLingShan Town, DongGuan City, Guangdong, China. Company Name: Wellkang Ltd Company Address: Suite B, 29 Harley Street LONDON W1G 9QR, England, United Kingdom Tel: +44 (20)30869438, 32876300 Fax: +44(20)76811874 Web: www.CEmark.com, www.CE-marking.com, www.CE-marking.eu, Email: AuthRep@CE-marking.eu

Company name: Babyono, Ul. Kowalewicka 13, 60-002 Poznań, Poland

#### EQUIPMENT SYMBOLS

|   |   |   |   |
|---|---|---|---|
|  | Refer to operating instructions   |  | Compliance with WEEE Standard   |
|  | This device complies with Part 15 of FCC( Federal Communications Commission) Rules.   |  | Indicates this device is in compliance with MDD 93/42/EEC. 0598 is the Notified Body Number |
|  | Manufacturer  |  | Restriction of Hazardous Substances   |
|  | Recyclable  |  | Authorized Representative in the European community   |
| <b>IP22</b>   | IP22: The first number 2: Protected against solid foreign objects ofØ12.5 mm and greater. The second number: Protected against vertically falling water |  | Serial number   |
|   |   |  | Warning   |

#### EMC DECLARATION

- 1) This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS; This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- 2) \* Caution: Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) \* Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4) \* Caution: this machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

#### GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC EMISSION

The Non Contact Infrared Body Thermometer is intended for use in the electromagnetic environment specified below. The customer of the user of the Non Contact Infrared Body Thermometer should assure that it is used in such an environment.

| EMISSION TEST   | COMPLIANCE     | ELECTROMAGNETIC ENVIRONMENT – GUIDANCE   |
|---|----------------|--|
| RF emissions<br>CISPR 11                                    | Group 1        | The Non Contact Infrared Body Thermometer use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                            |
| RF emission<br>CISPR 11                                     | Class B        | The Non Contact Infrared Body Thermometer is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions<br>IEC 61000-3-2                         | Not applicable |  |
| Voltage fluctuations/<br>flicker emissions<br>IEC 61000-3-3 | Not applicable |  |

#### GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Non Contact Infrared Body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of Non Contact Infrared Body Thermometer should assure that it is used in such an environment.

| IMMUNITY TEST   | IEC 60601 TEST LEVEL  | COMPLIANCE LEVEL<br>±6 kV contact                                | ELECTROMAGNETIC ENVIRONMENT<br>- GUIDANCE  |
|---|---|--|--|
| Electrostatic discharge<br>(ESD)<br>IEC 61000-4-2   | ±6 kV contact<br>±15 kV air   | ±15 kV air<br>Not applicable<br>Not applicable<br>Not applicable | Floors should be wood, concrete or ceramic tile.<br>If floor are covered with synthetic material, the relative humidity should be at least 30%.  |
| Electrical fast transient/<br>burst<br>IEC 61000-4-4  | ±2 kV for power supply lines<br>±1 kV for input/output lines  | 30A/m  | Mains power quality should be that of a typical commercial or hospital environment.  |
| Surge<br>IEC 61000-4-5  | ± 1 kV line(s) to line(s)<br>± 2 kV line(s) to earth  |  | Mains power quality should be that of a typical commercial or hospital environment.  |
| Voltage dips, short<br>interruptions and voltage<br>variations on power<br>supply input lines<br>IEC 61000-4-11 | <5% UT (>95% dip in UT)<br>for 0.5 cycle<br>40% UT (60% dip in UT)<br>for 5 cycles<br>70% UT (30% dip in UT)<br>for 25 cycles<br><5% UT (>95% dip in UT)<br>for 5 sec |  | Mains power quality should be that of a typical commercial or hospital environment. If the user of the Non Contact Infrared Body Thermometer requires continued operation during power mains interruptions, it is recommended that the Non Contact Infrared Body Thermometer be powered from an uninterruptible power supply or a battery. |
| Power frequency<br>(50Hz/60Hz) magnetic<br>field IEC 61000-4-8  | 30A/m   |  | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.  |

NOTE UT is the a.c. mains voltage prior to application of the test level.

## GUIDANCE AND MANUFACTURE'S DECLARATION – ELECTROMAGNETIC IMMUNITY

The Non Contact Infrared Body Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the Non Contact Infrared Body Thermometer should assure that it is used in such an environment.

| IMMUNITY TEST                 | IEC 60601 TEST LEVEL        | COMPLIANCE LEVEL | ELECTROMAGNETIC ENVIRONMENT - GUIDANCE  |
|-------------------------------|-----------------------------|------------------|---|
| Conducted RF<br>IEC 61000-4-6 | 3 Vrms<br>150 kHz to 80 MHz | Not applicable   | Portable and mobile RF communications equipment should be used no closer to any part of the Non Contact Infrared Body Thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  |
| Radiated RF<br>IEC 61000-4-3  | 10V/m<br>80 MHz to 2.5 GHz  | 10 V/m           | <p><b>RECOMMENDED SEPARATION DISTANCE</b></p> $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2,3\sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Non Contact Infrared Body Thermometer is used exceeds the applicable RF compliance level above, the Non Contact Infrared Body Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Non Contact Infrared Body Thermometer.

b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

**RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE NON CONTACT INFRARED BODY THERMOMETER.**

The Non Contact Infrared Body Thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Non Contact Infrared Body Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Non Contact Infrared Body Thermometer as recommended below, according to the maximum output power of the communications equipment.

| RATED MAXIMUM OUTPUT POWER OF TRANSMITTER (W) | SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER (M) |  |   |
|---|---|--|---|
|   | 150 KHZ TO 80 MHZ<br>$d = 1,2\sqrt{P}$                        | 80 MHZ TO 800 MHZ<br>$d = 1,2\sqrt{P}$ | 800 MHZ TO 2.5 GHZ<br>$d = 2,3\sqrt{P}$ |
| 0,01  | 0,12  | 0,12                                   | 0,23                                    |
| 0,1   | 0,38  | 0,38                                   | 0,73                                    |
| 1   | 1,2   | 1,2                                    | 2,3                                     |
| 10  | 3,8   | 3,8                                    | 7,3                                     |
| 100   | 12  | 12                                     | 23                                      |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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