# BEAUTIFLY

Prestige Tri-Action

Golden  
series

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Charger:	3.7V
Entry:	AC100-240V 50/60Hz
Exit:	DC5V 1A
Battery:	Built-in lithium-ion battery
timer:	15 minutes
Working time:	7-11 hours
Included:	Device, 3 attachments, charging station, power adapter, user manual, warranty card

### Security

Before using the product, read this manual and follow the instructions contained therein when using the device.

### Accessories [ see page 2, f1 ]

1. Brush attachment. Suitable for the care of the scalp, face and body. It uses electrode needles (32 pieces), LEDs
2. Face attachment. Firms and smoothes the skin of the face. It uses electrodes and LEDs
3. Body attachment. Firms and improves the appearance and condition of the skin. It uses an electrode section

### Usage Tips

- The device uses positive and negative electrodes to form a circuit that generates an electric current when in contact with the skin. For best results, make sure the electrodes are in full contact with the skin during the treatment.
- The sensation caused by electrical stimulation may be different on different parts of the body. Therefore, start the treatment from level 1 (when only LED 1 is lit) and adjust the power to ensure comfort of use.
- Do not use the device unreasonably or on skin areas with problems or abnormalities.
- Temporary redness may appear on the skin after using the device. If these symptoms do not improve over time, consult your doctor.
- Make sure the device is fully charged before use.
- Before using it for the first time or when the battery is low, charge it using the charging station.
- Before placing the charging terminal in the charging station, carefully dry it with a cotton swab and wait until it is completely dry.
- Full charge time is approximately 2.5 hours and operating time is approximately 7 hours when using the device at level 5 (with a brush attachment, at a room temperature of 27 degrees Celsius). However, actual operating and charging times may vary depending on previous battery level and charging conditions.

Before using the device, it is recommended to perform a patch test on the skin:

1. Apply an appropriate amount of EMS conductive gel or other cosmetic to clean skin.
2. Turn on the device and check its operation on the inside of the arm for about 1 minute using the lowest power level (when only LED 1 is lit).
3. If you experience any unusual reactions during the test, such as rash, itching or pain, stop the procedure immediately.

### Modes and functions:

The device uses two types of electrical stimulation that automatically adapt to the attachment you are using.

**Skin lifting mode\*** - through pulses sent radially from the electrode diodes, it stimulates the skin and the hair follicle muscles, providing lifting and a clear improvement in tension. Used with an overlay in the form of a brush for the care of the scalp, face, shoulders and neckline. Used dry, however, when using a moderate level of stimulation, it is possible to use it on damp skin or hair. It does not require prior application of gel or other EMS-conducting cosmetic.

Using the lifting mode on the scalp, face and body:

- Perform the treatment once a day for a maximum of 30 minutes, assuming that the session lasts from 5 to 15 minutes.
- Start with the lowest operating mode (level 1) to get used to the tingling sensation caused by the electrical stimulation.
- Do not use the device on the nose, eyelids, ears, mouth and eyes.
- For the best results, match the movements and the level of electrostimulation to the treated area of the body.

**Head care**

- Use on clean, combed hair and scalp
- In a slow zig-zag motion, comb your hair from the roots towards the back of the head
- From the back of the head, in a zig-zag motion, gradually move the device from above to the opposite clavicle
- Brush perpendicular to the skin, in full contact with the hairy part of the head; do not press the cap too hard to eliminate the risk of damage

**Face care:**

- Remove make-up and clean your face thoroughly
- Starting from the chin, lower jaw or nasolabial folds, gently touch the skin with the brush, moving towards the ears; perform similar movements on the forehead, starting from its base and moving towards the hairline; repeat the series 2 to 4 times

**Neck, shoulder and cleavage care:**

- Cleanse the skin thoroughly
- Move the brush in a smooth motion, starting from the ear area towards the neckline; from the back of the ear, go through the head towards the collarbone; Finally, move the device around your neck and shoulders.

\*During the procedure, you feel a tingling sensation caused by electrostimulation.

**Deep EMS mode\*** - through a deep EMS pulse, it strengthens the muscles and increases the production of collagen and elastin. Visibly reduces fine lines and wrinkles and improves skin tone. Stimulates blood and lymph circulation, giving the skin a young and radiant look. Used with overlays for the face and body, which ensure comfort of use. Requires the use of a gel or other cosmetic that conducts EMS currents. If the skin is dry, the electrostimulation may be difficult to feel. It is recommended to avoid the use of water or lotion as they can dry the skin and reduce the effectiveness of the electrical stimulation while causing discomfort.

**Using the EMS mode on face and body skin:**

- Perform the treatment once a day for a maximum of 30 minutes, assuming that the session lasts from 5 to 15 minutes.
- Start with the lowest operating mode (level 1) to get used to the tingling sensation caused by the electrical stimulation.
- Do not use the device on the nose, eyelids, ears, mouth, eyes and throat.
- To ensure comfortable use, do not press the device too hard against the cheeks, chin or other parts of the face when using it.
- Remember that when using the device, you may experience a stronger sensation of stimulation on the muscles of your head.

**face care**

- Before use, remove make-up and clean the skin thoroughly
- Use the facial care cap and apply gel or other EMS conductive cosmetic
- Starting from the chin, lower jaw or nasolabial folds, locally hold the device for 3-5 seconds, moving towards the ears; perform similar movements on the forehead, starting from its base and moving towards the hairline; repeat the series 2 to 4 times

**Neck and cleavage care**

- Thoroughly clean the skin and apply a gel or other EMS-conductive cosmetic
- Apply spotwise, for 3-5 seconds, holding the device in the selected place and moving it to other parts of the neck and cleavage; repeat the series 2 to 4 times

**body care:**

- Thoroughly clean the skin and apply a gel or other EMS-conductive cosmetic
- Note that it may be difficult to feel the EMS stimulation in this mode if the skin is dry
- Use on the muscles of the arms and legs by moving the device from bottom to top
- Apply on the belly in a zigzag motion, starting from the sides and moving towards the belly button
- Apply on the buttocks by moving the device from bottom to top
- Gently hold the device in the places you want to focus on
- Avoid stopping the device on bones

\*During the procedure, a slight tingling sensation is felt.

For more effective care, different attachments can be used during treatment. The lifting effect through electrical stimulation will ensure the use of a brush type attachment for 15 minutes a day to improve the appearance and tension of the skin before applying makeup. The recommended usage time for each area is as follows: - Face, head, shoulders and neck (5 minutes).

To exercise the facial muscles and further improve the condition of the skin, it is recommended to repeat the treatment with the facial attachment 3-4 times a week for 20 minutes before bedtime. The recommended application time for each area is as follows:

- Face, head, shoulders and neck (10-15 minutes) - Scalp (5-10 minutes) - Face, head, shoulders and neck

It should be remembered that the recommended time of using the device on individual parts of the body should be adapted to individual needs and the condition of the skin of the body, face and head. Use common sense and if you experience any discomfort or unusual sensations, stop using the device immediately and consult a physician.

#### Storage

- Keep the device dry and clean, out of the reach of children.
- Keep away from high temperature.
- Store the device in a cool, dry place.

#### Cleaning

- Use neutral detergent and a soft cloth for daily maintenance. Do not clean the device with organic solvents or other chemicals. Do not rinse with water.
- Avoid high-pressure washing and cleaning with hard objects that may scratch its surface.
- Between treatments on different parts of the face and body, disinfect the tips using agents intended for cosmetic devices

#### Landing:

- Connect to USB for charging. The device cannot be used while charging.
- When the battery level is low, LED 1 flashes quickly 3 times every 3 seconds.

#### Precautions and safety

- Use the device in a well-ventilated room with a temperature of 5°C~40°C and a humidity of 30%~85%.
- Keep the device away from electrical appliances, sources of high-frequency electrical radiation, etc., to avoid interference during use.
- Do not disassemble, repair or modify the device.
- When using for the first time, control the time (5-15 min) and pay attention to the frequency.

Excessive use may cause discomfort. If you feel uncomfortable, stop using the device.

- If the device does not work properly, stop using it immediately and contact the service center.

#### Warnings and contraindications

- If the device has not been used for a long time, wipe it down and turn it on to check that it works properly. Then charge the device. If any irregularities occur, please contact the service.
- Persons with implanted electronic devices, such as pacemakers, must not use the device.
- Using the device near a short wave or microwave processing device (eg 1M) may cause the output of the device to become unstable.
- Do not use the device during pregnancy and menstruation.
- Keep the device and its accessories out of the reach of children.
- Do not immerse the device in any liquid and do not wipe the product with a wet cloth.
- Turn off the device before charging.
- Improper use or adjustment of the device may void the service warranty.
- It is forbidden to use the device by people with wounds, fractures, scars, cancer, varicose veins and local poor skin elasticity.
- It is forbidden to use the device for people with serious heart, liver, kidney diseases or high fever.
- \*Scalp treatment should not be performed by people with hair extensions, people who have had hair transplants and people wearing certain types of hairpieces (such as adhesive tape or extensions).

More information about the product is available on the manufacturer's website at:

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EN: Electrical and electronic devices, as well as batteries and accumulators, must not be mixed with other household waste. This is indicated by the crossed-out trash can symbol on the packaging. Used electrical and electronic devices, batteries, and accumulators must be taken to local selective waste collection points or to the retailer, and more information can be found in your municipality. Electrical and electronic equipment, batteries, and accumulators may contain hazardous substances, mixtures, and components harmful to the environment and human health. Improper handling or damage to them can lead to health hazards or environmental pollution during further disposal or recycling.

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EN: The manufacturer declares under his sole responsibility that the product Prestige Tri-Action Golden is in conformity with the European Parliament and Council Directive 2014/30/EU of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility, the European Parliament and Council Directive 2014/35/EU of 26 February 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits, and the European Parliament and Council Directive 2011/65/EU of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment, together with the Commission Delegated Directive (EU) 2015/863 of 31 March 2015 amending Annex II to the European Parliament and Council Directive 2011/65/EU with regard to the list of restricted substances. The full text of the declaration of conformity is available at the following internet address:

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