USER MANUAL

* Images and diagrams referenced in this manual are for illustrative purposes only, The actual appearance may differ.

The Instructions can be found at www.skg.com

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
EN

Heat level

Heat level	- Low heat - Medium heat - High heat
Attention	To avoid low-temperature burns, please carefully use the heat level. If you experience any discomfort, please take a break promptly.

Massage speed

Low speed	Slow rotation massage with the massage head
High speed	Fast rotation massage with the massage head

Massage mode

Kneading mode	For individuals with discomfort in the neck, shoulders, and back without applying too much pressure, this mode uses low-speed, gentle kneading to ease discomfort, with a 10-minute cycle per session.
Pinching mode	This mode employs a rhythmic, progressive massage technique to enrich the user experience and alleviate muscle tightness, with a 10-minute cycle per session.
Acupressure mode	Utilizing strong tuina (Chinese therapeutic massage) techniques, this mode addresses issues of soreness, swelling, and pain, with a 10-minute cycle per session.

Charging/Power Indicator Light Instructions

EN

Light Color	Light Status	Description
Green Light	Steady On	Battery capacity > 25%
Orange Light	Steady On	Battery capacity < 25%
Orange Light	Breathing	Battery charging
Green Light	Breathing	Battery fully charged

Troubleshooting

EN

Status	Possible Cause/Handling Method
Cannot Start	- Battery is low, charge the product before use Ensure the main power button on the product is pressed
Cannot Charge or No Charging Light Displayed	- Check if the charging cable is properly connected and if the charger specifications meet the standard requirements

Product Wearing Instructions

EN

1. Please adjust the strap first. The tightness of the strap affects the massage intensity. Please adjust the strap to the appropriate length (ensure that the length of the straps on both sides are consistent).
2. Place the machine on your neck.
3. Please pass the straps on both sides from the front of the body under the armpits. Then pull the buckle to the back and insert it.

After use

EN

1. After use, please press the buckle to release the shoulder strap.
2. Loosen the back buckle.
3. Please carefully remove the product from the shoulders to prevent dropping.

Notes for Product Scrap

Please support the recycling of rechargeable battery (lithium-ion battery).

 Danger	The battery in this product is a dedicated rechargeable battery, which is not intended for other products. Please do not charge the removed rechargeable battery. - Do not throw it in the fire, or heat it up to high temperature. - Do not pierce it with nails, impact it by gravity or decompose it for modification. - Do not connect the positive and negative poles of the battery directly with metal. - Do not store it together with necklaces, hairpins, etc. - Do not charge, use or put in the places in high temperature such as fireside or in the blazingsun. - The situations above may cause heating, fire and breakage.
 Warning	Keep the removed rechargeable battery out of the reach of infants and young children in case that it is swallowed by mistake. When the device is mistakenly swallowed, please consult a doctor immediately. When the battery fluid leaks, do not touch the fluid by hand, and it can be processed in the following way. Do not rub your eyes with your hands when battery electrolyte enters your eyes in case of blindness. Rinse your eyes with clean water immediately and seek medical advice. When your body or clothing is contaminated by the battery electrolyte, it may cause skin irritation or injury. Please rinse with water. If you still feel discomfort, seek medical advice immediately. Only the scrapped products can be decomposed. Otherwise, the accident of fire, electric shock, injuries may be resulted in. When you need to repair the battery (not to replace the battery), please consult the dealer where you purchased this product.

The lithium-ion battery in this product is a valuable resource that can be recycled and reused, which is not allowed to be replaced. Please hand it over to a professional agency that can recycle the rechargeable batteries when this product is scrapped. Before scrapping this product, remove the battery from this appliance. Before removing the battery, the appliance must be powered off. The battery should be disposed of safely, which can be taken out in the steps below:

1. With professional tools, open the upper and lower cases.
2. Remove the circuit board and the motor.
3. Take out the battery, cut the battery wire with professional pliers, and then take out the battery.

Steps for removal of battery:

1. Open the shell with tools;
2. Loosen the screws on the circuit board;
3. Cut the battery lead off with scissors;
4. Take out the battery;

Precautions

Usage Restrictions and Locations

Please do not use this product in conjunction with the following medical electronic devices to avoid interference and potential harm:

- Internal medical electronic devices such as pacemakers
 - Life-sustaining medical electronic devices such as artificial heart-lung machines
 - Medical electronic devices such as electrocardiographs
- Do not use it in flammable or explosive environments such as petrol stations, high dust workshops, or gas stations.

The appliance has a heated surface. Persons insensitive to heat must be careful when using the appliance. This appliance contains batteries that are non-replaceable.

This product is not a medical device and does not have the function to diagnose or treat neck diseases. Patients with neck conditions should consult medical personnel to avoid delaying treatment. Do not use this product while driving, climbing, or performing other hazardous activities, as it may cause accidents.

It is recommended that the product be used for no more than 20 minutes (2 cycles) at a time to avoid the risk of low-temperature burns and discomfort due to excessive massage.

Start with a low speed and gradually increase to a higher setting according to personal comfort to avoid discomfort or injury from high-intensity settings.

Do not use it while sleeping or in a state of impaired consciousness.

Individuals with cerebrovascular diseases or other circulatory disorders, inflammation, swelling, suppurative in the neck area, or cervical vertebrae diseases should not use this product, as it may exacerbate symptoms.

Keep the product out of reach of infants, children, and anyone who is fully or partially unconscious. Do not allow them to use the product to prevent accidents.

Caution is advised for users over the age of 70 (use only with a doctor's permission and guidance). For users over the age of 65, due to more delicate skin and less sensitive nerve endings, please use a lower setting and strictly monitor usage time.

Warnings

This product is a personal care massage device and is not a substitute for medication or medical devices. Not for use by children or pregnant women.

Women during menstruation should use with caution based on their physical condition.

Contains a lithium battery; do not place in high-temperature environments.

Individuals with sensitive skin, facial dermatitis, or facial implants should use under medical guidance.

Users with skin allergies to metals or heat should exercise caution when using this product.

This product has a heating function; users insensitive to heat should use the low heat setting and strictly control usage time.

Discontinue use immediately if discomfort occurs during use.

Do not attempt to disassemble the product.

Do not attempt to replace the battery yourself.

Individuals with cervical spine injuries should consult healthcare professionals before using this product.

The device has a heating surface; users insensitive to heat must be cautious.

Use in a normal environment: Temperature: 0°C-40°C; Relative Humidity: <80% RH, non-condensing; Atmospheric Pressure: 50Kpa-106Kpa.

Choose normal transportation and storage conditions: Temperature range: -20°C to +55°C; Relative Humidity range: <95% RH, non-condensing; Atmospheric Pressure: 50Kpa-106Kpa.

This product contains a lithium battery; do not disassemble, modify, immerse in water, or throw into fire as it may cause damage or explosion.

Do not apply pressure to the heart during use, as it may cause accidents or discomfort.

Do not use with medical devices or more than one of this product simultaneously, as it may cause nausea or discomfort.

Do not use for purposes other than those specified in the instruction manual, as it may cause accidents.

Do not disassemble, repair, or modify the product yourself, as it may lead to malfunction or electric shock.

For an optimal experience, do not forcefully pull or lean on the product!

Consult a doctor before use in the following situations, as there may be a risk of accidents or discomfort:

Patients currently under medical treatment or feeling unwell, including those with infectious diseases

Patients with heart or neurological abnormalities

Patients with sensory disturbances or abnormalities in the skin

Patients with a body temperature exceeding 38 degrees Celsius (feverish)

Patients currently taking medication

Patients with abnormal blood pressure

Patients unable to express their own will

Patients with sensory disturbances due to severe peripheral circulatory disorders caused by conditions such as diabetes

Do not use in the following situations:

- Patients with malignant tumors
- Near the heart area of the human body
- Patients with metal implants
- Pregnant or postpartum women
- Users of pacemakers
- Patients with acute suppurative inflammation
- Patients with coagulation disorders
- Patients with a tendency to bleed
- Patients with spastic paralysis
- Patients with severe cognitive impairments
- Patients in the spastic phase following cerebrovascular accidents
- Individuals with abnormal neurogenic stress responses
- Individuals with unclear consciousness or limb joint contracture deformities

Precautions:


- Please read the instruction manual thoroughly before using this product.
- If you feel unwell due to an abnormality with the product, stop using it immediately. If you notice any unusual sensations in your body or skin, consult a doctor and follow their advice.
- In the event of a malfunction that prevents operation, immediately disconnect the power and remove the product to prevent the risk of overheating or short-circuiting. Contact customer service for assistance.
- Do not plug or unplug the charging cable with wet hands, as this may lead to electric shock or accidents.
- Do not wear necklaces while using the product.
- Do not use the product while lying flat.
- Wear the product correctly, positioning it on the cervical spine area and ensuring symmetrical placement on both sides to avoid accidents or discomfort.
- When disconnecting the charging cable from the product, do not pull on the cord; instead, grasp the plug to remove it to prevent electric shock or damage to the product.
- When disconnecting the charging cable from the product, do not pull on the cord; instead, grasp the plug to remove it to prevent electric shock or damage to the product.
- If using the product for the first time or after a long period of non-use, charge the battery according to the instructions to ensure it functions properly.
- It is normal for the product's surface to become warm during extended operation.
- The product will massage the neck area appropriately; do not massage continuously for extended periods. If discomfort occurs, take a break before using it again.

Suggestions:

- Ensure that the main unit buttons are functioning properly before use.
- To prevent pinching injuries, be cautious to avoid placing fingers or other body parts in narrow parts of the device.
- When disposing of this product, please do so in accordance with local environmental protection regulations to prevent environmental pollution.
- Recommended for use by individuals aged 16 and above.

Please do not use or store this product in the following locations:

- Places exposed to direct sunlight.
- Areas with high temperatures and humidity.
- Locations with the presence of water.
- Dusty environments.
- Near open flames.
- In strong electromagnetic fields.
- Places where the product may be subject to vibrations or shocks.

This device in compliance with the essential requirements and other relevant provisions of Directive (EU)2023/988.

CAUTION : RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.

DISPOSE OF USED BATTERIES

ACCORDING TO THE INSTRUCTIONS

EU Compliance Statement:

SKG Health technology Co.,Ltd. Here by declares that this device is in compliance with the essential requirements and other relevant provisions of the Directive (EU)2023/988. A copy of the EU Declaration of Conformity is available online at <https://www.skg.com/pages/skg-declaration-of-conformity>.

Old electrical appliances must not be disposed of together with the residual waste, but have to be disposed of separately. The disposal at the communal collecting point via private persons is for free. The owner of old appliances is responsible to bring the appliances to these collecting points or to similar collection points. With this little personal effort, you contribute to recycle valuable raw materials and the treatment of toxic substances. Caution:

- If you experience skin irritation when using the body massager, please refrain from using it, and wait two to three days to see if your symptoms ease. If symptoms persist or worsen, please consult a physician.
- If you have eczema, allergies or asthma, you may be more likely to experience a skin irritation or allergy from using the massager.
- The appliance it must only be supplied at safety extra low voltage corresponding to the marking on the appliance.
- This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved.
- Children shall not play with the appliance.
- Cleaning and user maintenance shall not be made by children without supervision.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Technical parameters


Product name	Neck & Shoulder Massager	Rated input power	10W
Model	H7 Ultra E	Rated input voltage	5V=
Battery capacity	2500mAh/18.5Wh		

Note: This product is rechargeable and is suitable for charging with typeinterface. Users need to charge with a standard charger with an output of 5V = 2A

SUCCESS COURIER SL
CALLE RIO TORMES NUM.1, PLANTA 1, DERECHA,
OFICINA 3, Fuenlabrada, Madrid, 28947 Spain
succeservice2@hotmail.com

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SKG Health Technologies Co., Ltd.